

## Potentiometric Error Detector



**Product Categories:** [Control Lab](#), [Control Lab](#), [Electronics](#), [Engineering Equipment](#)

**Product Page:** <https://www.labappara.com/product/potentiometric-error-detector/>

### Product Description

Potentiometric Error Detector

#### **TECHNICAL SPECIFICATIONS**

To study of potentiometer as an error detector.

Two Full 360° rotation servo-potentiometers with 1° degree resolution.

IC regulated DC excitation for both potentiometers.

AC excitation at 400 Hz approx.

Built in demodulator circuit.

3½ digit DVM for carrying measurements.

IC regulated power supply for rest of the circuitry.

220V ±10%, 50Hz mains operated.

User's manual and patch cords.