

Ganong's Potometer



Product Categories: [Plant Physiology Apparatus](#)

Product Tags: [Laboratory Equipment](#)

Product Page: <https://www.labappara.com/product/potometer-ganong/>

Product Description

Ganong's Potometer - Laboratory Equipment

A Ganong's potometer is a piece of apparatus used to measure the rate of water loss from a plant.

The potometer measures the water lost through transpiration of the plant and not the water taken up by the plant.

Borosilicate glass for measuring absorption of water by transpiring plants.

It is used for measuring the rate of water uptake by a leafy plant shoot.

Ganong's Potometer that are fabricated, using exceptional quality imported glass tubing that offers durable structure to them.

Features:

High efficiency

Sturdy construction

Optimal performance

Ganong's Potometer on stand with packed in card box.

Ganong's Potometer with borosile glass.