

## Projection Microscope 6"



**Product Categories:** [Microscope](#), [Microscopes](#)

**Product Tags:** [biology lab equipment](#)

**Product Page:** <https://www.labappara.com/product/projection-microscope-6/>

### Product Description

#### Projection Microscope 6" - Scientific Equipment

Deluxe with two objective and one eyepiece, fitted with 12 V-21 W halogen lamp with intensity control knob. The structure of microscopic specimens becomes visible to the naked eye by taking the projection on the screen, during observation and demonstration to a group of students. The microscopes are compact in design and can easily contrast and distortion free magnified images on a fine ground glass screen. Well packed in plywood box.

Highly sensitive coarse and fine motion, triple revolving nose piece with 2 objectives 10 x, 45 x & 10 x eye piece with 6" Projection Dome in Thermocol box.

#### **Features:**

Long working life

Low maintenance

Hassle free functionality

Simple operation

High strength

Fine finish