


## Quartz Cell with PTFE Cover



**Product Categories:** [Laboratory Glassware](#), [Quartz Ware](#)

**Product Tags:** [glassware](#)

**Product Page:** <https://www.labappara.com/product/quartz-cell-ptfe-cover/>

### Product Description

Quartz Cell with PTFE Cover - Laboratory Glassware Equipment

All cuvettes are made with modern precision manufacturing methods, it is no longer necessary to pre-match cuvettes.

#### **Features:**

Highly accurate results, only need to test sample one time

Accurate repeat-ability, Precision Cells cuvettes are made with only the highest quality quartz and glass materials

Hand made, ensuring the highest quality optical polish with no defects.

Every cell is thoroughly inspected before leaving the factory.

Precision Cells cuvettes are guaranteed to give you the best optical transmission.

Fire fused to ensure the best bonding.

Three premium materials available: (Optical Glass, Pyrex Glass (10mm only), UV Quartz, IR Quartz)

Lightpaths from 1 mm up to 100 mm available.

Two windows optically polished.

Includes PTFE cover.

Place your order on-line using the chart below.

**Recommended accessories:**

Spacer for adapting 1, 2 and 5 mm versions to a standard 10 mm cuvette holder:

A25 Cuvette Spacer

Spring-based mount for fitting any versions less than 5 mm to a standard 10 mm cuvette holder: A2o Series: Cuvette Mount

Cuvette Holder for cells with 1-100 mm path-length: CHR1/x

Mixer for use with 10 mm cells: P69 Mixer

Magnetic stirrer for use with 10 mm cells: P73 Stirrer

Cleaning swabs: SwabMFC

Spare PTFE covers

10 mm path, 3.5 ml in pair.