


## Rack For Microcentrifuge Tube



**Product Categories:** [Centrifuges Equipment](#), [Laboratory Equipments](#)

**Product Tags:** [biology lab equipment](#), [chemistry lab instruments](#), [physics lab instruments](#)

**Product Page:** <https://www.labappara.com/product/rack-micro-centrifuge-tube/>

### Product Description

Rack For Microcentrifuge Tube - Scientific Laboratory Equipment

Moulded in Poly-carbonate, these autoclavable racks are made for holding 1.5 ml. micro centrifuge tubes. These racks have a unique stacking feature which prevents movement in any direction. The holes are numbered & lettered for easy identification.

Each hole has molded-in-numbers for easy sample identification.

Exclusive and best range

Very well tested

Efficient

### Features

Sturdy rack with rubbered support at the bottom for stability

Racks lock together on each side to allow user to attach several for larger experiments

Compatible with both 1.5 ml and 2.0 ml Micro Centrifuge Tube

Easy removal and insertion of tubes.

Each hole has molded-in-numbers for easy sample identification

Can be autoclaved.