

## Rack Scintillation Vial



**Product Categories:** [Laboratory Plastic Ware](#), [Vails](#)

**Product Tags:** [biology lab equipment](#), [chemistry lab instruments](#)

**Product Page:** <https://www.labappara.com/product/rack-scintillation-vial/>

### Product Description

#### Rack Scintillation Vial - Scientific Laboratory Equipment

These racks for scintillation vials, moulded in polypropylene, are autoclavable & corrosion resistant. These is an individual rack for both 8 ml & 20 ml vials. These rugged & stable racks can hold scintillation vials for convenient storage, handling & sample preparations.

**Material :** Polypropylene

### **PARTICULARS PACKING (Pcs)**

Rack for 90 vials of 8ml capacity

2

Rack for 50 vials of 20ml capacity

2

Product data have been exported from - [Labappara scientific instrument](#)

Export date: Wed Nov 27 3:32:56 2024 / +0000 GMT