

Rain Gauge



Product Categories: [Measurements](#), [Physics Equipments](#)

Product Tags: [physics Equipments](#)

Product Page: <https://www.labappara.com/product/rain-gauge/>

Product Description

Rain Gauge - Physics Laboratory Equipment

Constructed from copper. It has funnel, an inner receiving vessel and an outer body all of copper, and a glass measuring cylinder. The set of rain gauge is provided with metal case having removeable lid to collect rainfall and a plastic calibrated rain gauge jar duly marked in centimeters.

Features:

Quality tested

Hassle free performance

High strength