


## Reagent Reservoir



**Product Categories:** [Laboratory Plastic Ware](#), [Storage Boxes Carriers](#)

**Product Tags:** [biology lab equipment](#), [chemistry lab instruments](#), [Manufacturer of Laboratory equipment](#), [math lab instruments](#), [physics lab instruments](#), [science lab equipment](#)

**Product Page:** <https://www.labappara.com/product/reagent-reservoir/>

### Product Description

Reagent Reservoir - Scientific laboratory equipment

**Reagent Reservoirs:-** . Made of polypropylene these reagent reservoirs are repeatedly autoclavable and can be used with all multi-channel pipettes. The Universal Reagent Reservoir is a graduated single reagent reservoir with maximum volume of 50 ml, reservoir with 8 separate numbered 5 ml. Reagent channels & a reservoir with 12 separate numbered 5 ml. Moulded in Polypropylene, this autoclavable reagent reservoir is helpful for multi-channel pipettes to be used into micro titer plates.

**Material :** Polypropylene

PARTICULARS      PACKING (Pcs)

50 ml                      6 Pcs

75 ml                      12 Pcs