


Reagent Reservoir



Product Categories: [Laboratory Plastic Ware](#), [Storage Boxes Carriers](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#), [Manufacturer of Laboratory equipment](#), [math lab instruments](#), [physics lab instruments](#), [science lab equipment](#)

Product Page: <https://www.labappara.com/product/reagent-reservoir/>

Product Description

Reagent Reservoir - Scientific laboratory equipment

Reagent Reservoirs:- . Made of polypropylene these reagent reservoirs are repeatedly autoclavable and can be used with all multi-channel pipettes. The Universal Reagent Reservoir is a graduated single reagent reservoir with maximum volume of 50 ml, reservoir with 8 separate numbered 5 ml. Reagent channels & a reservoir with 12 separate numbered 5 ml. Moulded in Polypropylene, this autoclavable reagent reservoir is helpful for multi-channel pipettes to be used into micro titer plates.

Material : Polypropylene

PARTICULARS	PACKING (Pcs)
50 ml	6 Pcs
75 ml	12 Pcs