

Reciprocity Theorem



Product Categories: [Electronics](#), [Engineering Equipment](#), [Network Theorem's Verifications \(In DC Circuits\)](#), [Network Theorem's Verifications \(In DC Circuits\)](#)

Product Page: <https://www.labappara.com/product/reciprocity-theorem/>

Product Description

Reciprocity Theorem

Analog Meters:

Ammeter 30mA DC.

Voltmeter 15V DC.

Power Supplies:

DC Supply IC Regulated 0-15V DC, 150mA.

Operated on Mains power 230V, 50Hz +10%

Components are mounted on the panels are:

Resistances 2Watts

Variable Resistor

Voltage Control through Potentiometer.

FEATURES:

Front panel built with high class insulated Printed Circuit Board sheet with well printed circuits and symbols.

Fuse for Short Circuit protection

Instruction manual.

Connections are brought out through 4mm Colored Sockets.

Patch Cords 4mm.

The trainer is housed in ABS Plastic cabinet.

Size of the trainer set 12"x8"