

## Reduction Adapters



**Product Categories:** [Adapters Borosilicate Glass](#), [Laboratory Glassware](#)

**Product Tags:** [glassware](#)

**Product Page:** <https://www.labappara.com/product/reduction-adapters/>

### Product Description

Reduction Adapters - Laboratory Glassware Equipment

All reduction adapters are made from borosilicate glass and joints are interchangeable with other apparatus. An outer interchangeable joint at the top (Socket) and a smaller inner interchangeable joint (Cone) at the base.

Available in socket size - 14/23, 19/26, 24/29, 29/32, 34/35, 40/38

**Features:**

Chemical resistance

Leakage proof

Flawless finish