

REGNAULT'S CALORIMETER



Product Categories: [Heat And Thermodynamics](#), [Physics Equipments](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/regnaulsts-calorimeter/>

Product Description

Regnault's Calorimeter - Physics Apparatus

Renault's calorimeter is an elegant apparatus to secure accurate results in experiments on heat.

This is first described in 1840.

With it one can measure the amount of heat energy required to change the temperature of a solid body by a known amount.

This knowledge has both practical and conceptual value.