

## Resistance Box Dial Type



**Product Categories:** [General Physics](#)

**Product Tags:** [electronics](#), [physics lab instruments](#)

**Product Page:** <https://www.labappara.com/product/resistance-box-dial-type/>

### Product Description

Resistance Box Dial Type - Physics Laboratory Equipment

Designed for quick and accurate work.

The resistor values shown on thick Bakelite sheet and read by a pointer.

Resistance boxes each dial consists of 10 steps.

Constant coils accuracy  $\pm 0.5\%$  to  $\pm 0.1\%$

And with manganin coil accuracy  $\pm 0.2\%$  to  $\pm 0.05\%$

Comprises circular dials mounted underneath a metal cover plate, inside sturdy aluminium extrusion case.

Each dial has 10 resistances, with a knob at the top which rotates with positive click stop to select the desired value as indicated alongside index mark on the knob.

Provided with a pair of 4mm sockets for connecting the resistances across a circuit and a yellow socket for earthing.

Multiple dials are connected in series to add to the total resistances of each dial

Available in - One dial, two dial, three dial, four dial, five dial