


Resistance Box Plug Key



Product Categories: [Educational equipments](#), [General Physics](#)

Product Tags: [chemistry lab instruments](#), [Laboratory Equipment](#), [Scientific Equipment](#)

Product Page: <https://www.labappara.com/product/resistance-box-plug-key/>

Product Description

Resistance Box Plug Key:

Teak wood box fitted with 7/16" rectangular brass, blocks accuracy + 0.1%

Provided resistance box is designed and manufactured by our deft professionals using premium quality components and advanced technology at our well-equipped production unit.

Resistance box in different sizes, models, capacities and configurations.

Offered resistance box is provided with multiple resistors & switches to ensure appropriate flow of current, by adjusting the resistor value.

This resistance box is available in the market at the most reasonable price.

Features:

Hassle-free operations

Long service life

Anti-corrosive body

High durability

Sturdiness

Perfectly designed

Available In Range

1 to 10 or 1-1

1 to 50 or 1-5

1 to 100 or 1-10

1 to 500 or 1-50

1 to 1000 or 1-100

1 to 5000 or 1-500

1 to 10000 or 1-1000