

Resistance Coil Transparent



Product Categories: [General Physics](#)

Product Tags: [electronics](#), [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/resistance-coil-transparent/>

Product Description

Resistance Coil (Transparent) - Physics Laboratory Instrument

Resistance coil in transparent plastic tube.

Accuracy + 0.2 any value 0.1 to 100 ohms is available.

Resistance coil in round Bakelite case.

Resistance with brass screws terminals.

Resistance coil enclosed in Bakelite tough case.

Non-inductively wound with double silk covered constant an coils on plastic rod enclosed in transparent plastic tube with two terminal and values figured.

Available in value - 1 k, 2 k, 5 k, 10 k, 15 k (k stand for 1000 ohms).