

## Resistance Coil



**Product Categories:** [General Physics](#)

**Product Tags:** [Physics Laboratory Equipment](#)

**Product Page:** <https://www.labappara.com/product/resistance-coil/>

### Product Description

Resistance Coil - Physics Laboratory Equipment

Resistance Coils are commonly used in Physics, Electronics Laboratories

Resistance coil in round Bakelite case with 4mm socket captive head terminals.

Resistance with brass screws terminals.

And then any value 0.1 to 100 ohms is available.

These resistance coils have round case and captive head terminals.

These resistance coils are able to deal with the resistance and rust efficiently.

#### **Features:**

Finished

smoothly

Good quality