

## Resistance In Series & Parallel Apparatus



**Product Categories:** [Characteristics And Application](#), [Characteristics And Application](#), [Electronics](#), [Engineering Equipment](#)

**Product Page:**

<https://www.labappara.com/product/resistance-in-series-parallel-apparatus/>

### Product Description

Resistance In Series & Parallel Apparatus

#### **Features :**

Instrument comprises of Two moving coil meters to measure voltage & current, Circuit diagram for series & parallel connections Printed & Resistances mounted on the front panel.