

## Resonance Apparatus (Export Quality)



**Product Categories:** [Physics Equipments](#), [Sound And Waves Oscillations](#)

**Product Tags:** [physics lab instruments](#)

**Product Page:** <https://www.labappara.com/product/resonance-apparatus-export-quality/>

### Product Description

#### **Resonance apparatus:**

Resonance apparatus consists of water container, transparent tube, metal tube of different metal.

The transparent helps to measure the water level against the adjacent scale.

Adjusting the water container height the water level also gets adjusted accordingly.

#### **Features**

Steel C.P tube and brass reservoir on heavy iron base.

Export quality brass C.P brass tube and brass reservoir.

Good quality & superior