

Resonance Apparatus



Product Categories: [Physics Equipments](#), [Sound And Waves Oscillations](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/resonance-apparatus/>

Product Description

Resonance Apparatus:

Resonance apparatus consists of water container, transparent tube, metal tube of different metal.

The transparent helps to measure the water level against the adjacent scale.

Adjusting the water container height the water level also gets adjusted accordingly.

Features

Steel C.P tube and brass reservoir on heavy iron base.

With stainless steel tube and brass reservoir

Good quality & superior

It Is used in School, College, University, Research and Other Laboratories.