

# Respirometer



**Product Categories:** [Plant Physiology Apparatus](#)

**Product Tags:** [Laboratory Equipment](#)

**Product Page:** <https://www.labappara.com/product/respiraometer/>

## Product Description

Respirometer - Laboratory Equipment

A Respirometer is a device used to measure the rate of respiration of a living organism by measuring its rate of exchange of oxygen and/or carbon dioxide. These respirometers are widely used to demonstrate and measure the oxygen uptake in the respiring material such as germinating seeds. The respirometers are provided with a glass vessel graduated 75 x 100 x 1 mm. along with with levelling tube. These are available complete with a stand and clamp.

### Feature

Respirometer is complete with graduated and levelling tubes.

Respirometer is made consists of rubber tubing and stand.

And then with in borosile Glass.

Rigid structure

Dimensionally accurate

Long service life