

Retort Ring



Product Categories: [Chemistry Equipments](#), [Lab Supplies](#)

Product Tags: [chemistry lab instruments](#), [physics lab instruments](#), [science lab equipment](#)

Product Page: <https://www.labappara.com/product/retort-ring/>

Product Description

Retort Ring - Laboratory Equipment

Retort ring made of mild steel rod, powder coated. Attached with machine pressed MS boss head, able to hold 16mm dia rod. Boss-head provided with plated thumb metal screw having moulded plastic head for securely clamping the bosshead to retort stand rod. All metal parts zinc plated.

Available in size - 2", 3", 4"

Feature -

Reliability

Fine finish

Dimensional accuracy

Corrosion resistance

Precisely designed

Good Quality.

Used in laboratories