

Retort Stand



Product Categories: [Chemistry Equipments](#), [Lab Supplies](#)

Product Tags: [chemistry lab instruments](#), [science lab equipment](#)

Product Page: <https://www.labappara.com/product/retort-stand-2/>

Product Description

Retort Stand - Scientific Equipment

A Retort Stand sometimes called a ring stand, is a scientific equipment which clamps can be attached to hold [test tubes](#) and other equipment such as [burette](#) which are most often used in titration experiments. Retort stand is one of the general lab equipment or metal wares being used to help with holding onto the bench equipment.

Available in size -

A

6"X 4"X 18"

B

7"X 5"X 24"

C

8"X5"X 24"

D

9"X 6"X 30"

Feature -

Heavy cast Iron base.

Hammertone paint.

With Zinc plated rod.