

Retort Stand



Product Categories: [Clamps, Holders Racks Stands](#), [Laboratory Plastic Ware](#)

Product Tags: [Laboratory Equipments](#), [Laboratory plastic ware](#)

Product Page: <https://www.labappara.com/product/retort-stand/>

Product Description

Retort Stand - Scientific Laboratory Plastic ware

A retort stand, sometimes called a ring stand, to which clamps can be attached to hold test tubes and other equipment such as burettes which are most often used in titration experiments. Due to its non-corrosive nature it is a better replacement for the metal ones.

These polystyrene/Polypropylene non-corrosive, moulded retort stands provide excellent replacement for the metal ones. Even the metallic rod has a polypropylene covering making it non-corrosive. The rod measuring 12.5 mm x 750 mm can be screwed in to the heavy base giving way to a perfectly stable retort stand.

Particulars	Packing
225x150 mm-Side Hole	1 Pc
300 x 200 mm-Centre Hole	1 Pc
300x200 mm-Side Hole	1 Pc