

Rheostats Grooved Type



Product Categories: [General Physics](#)

Product Tags: [electronics](#), [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/rheostats-grooved-type/>

Product Description

Rheostats Grooved Type - Physics Laboratory Equipment

Grooved tube wound with eureka wires on 5.6 cm 2 ½" dia tube in grooved to avoid loosening of wire.

Brushes of hard laminated brass are phosphorus bronze strips.

Finished in tident grey and bright chrome.

Rheostats are available in several ranges.

Temperature resistance 1000c above room temperature.

These are made with special non breakable vitreous enameled steel tubes.

The winding is locked into place with ceramic cement.

The tube is carried upon robust enameled die-cast end stands and the heavy duty sliding contact is of the multi leaf phosphor bronze type.

Available in length -

LENGTH

IN (CM)

15cm

6"

20cm

8"

25cm

10"

30cm

12"

35cm

14"

40cm

16"

45cm

18"

50cm

20"