

Rheostats



Product Categories: [Electricity](#), [Measurement & Circuits](#), [Physics Equipments](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/rheostats/>

Product Description

Rheostats - Physics Laboratory Equipment.

Rheostats 4.3 cm (1 ¾") dia tube wound with heavily oxidized eureka wire, mounted feet

Rheostats with laminated brush contacts of phosphorus bronze strips brass c.p guide rod.

Lock type terminal , brass parts and metal feet.

Available in length - 15 cm, 20 cm, 25 cm, 30 cm, 35 cm, 40 cm, 45 cm, 50 cm.

Available in cm - 6" 8" 10" 12" 14" 16" 18" 20".