

Ria Vial Screw Cap



Product Categories: [Laboratory Plastic Ware](#), [Vials](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/ria-vial-screw-cap/>

Product Description

Ria Vial Screw Cap - Scientific Laboratory Equipment

These Polypropylene Vials are unique because the cap fitted with an O-ring can be screwed into the vial making it a leakproof vial which can be used in RIA, Coagulation & Bacteriology.

The polypropylene tubes are clear and can withstand 3000 RCF during centrifugation while the transparent polystyrene tubes can withstand 1400 RCF. The PP tubes can accept common acids, alkalies & solvents at room temperature while the PS tubes can accept mild bases, weak acids but not organic solvents.

Particulars

Packing

12 X 75 mm/ PS

500 Pcs

12 X 75 mm/ PP

500 Pcs

Cap for 12 mm RIA Tube

100 Pcs