

RIA Vial



Product Categories: [Laboratory Plastic Ware](#), [Vails](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/ria-vial/>

Product Description

RIA Vial-Scientific Laboratory Equipment

These Polypropylene/Polystyrene tubes are generally used in RIA, Coagulation & Bacteriology. The Polypropylene tubes are clear and can withstand 3000 RCF during Centrifugation while the transparent Polystyrene tubes can withstand 1400 RCF. The PP tubes can accept common acids, alkalies & solvents at room temperature while the PS tubes can accept mild bases, weak acids but not organic solvents.

Material : Polypropylene / Polystyrene

Particulars	Packing
12 x 75 mm / PS	500 Pcs.
12 x 75 mm / PP	500 Pcs.
Cap. for 12 mm RIA Tube	100 Pcs.