

Ripple Factor Apparatus (Half Wave, Full Wave & Bridge Rectifier)



Product Categories: [Characteristics And Application](#), [Characteristics And Application](#), [Electronics](#), [Engineering Equipment](#)

Product Page: <https://www.labappara.com/product/ripple-factor-apparatus-half-wave-full-wave-bridge-rectifier/>

Product Description

Ripple Factor Apparatus (Half Wave, Full Wave & Bridge Rectifier)

Features :

Instrument comprises of AC Power Supply, 2 meters to measure Output Voltage & Output current, 4 PN Junction Diodes, Filter Circuit Kit, load Resistances selectable using band switch.