

Rocking Microtome



Product Categories: [Biology Lab](#), [Educational equipments](#), [Microtome](#)

Product Tags: [biology lab equipment](#)

Product Page: <https://www.labappara.com/product/rocking-microtome/>

Product Description

Rocking Microtome -Biology Laboratory Equipment - Cambridge rocking microtome which has given outstanding service to numerous micro chemists - It is designed for sectioning comparatively soft, wax-embedded tissues with a thickness range. Calibrated in, the section cutting thickness can be set easily and quickly by a single control. - It is recommended for elementary biology & pathology students. - The rocking microtome is an extremely reliable instrument comprising of three moving parts. - The lead screw mechanism incorporates quality precision engineering. - The cutting stroke is controlled by a single control lever which minimizes the operator skill required. - It carries wide choice of specimen holders cylindrical holders, standard orientating & wood block holders for sectioning was embedded specimens.