


Rockwell & Brinel Hardness Tester



Product Categories: [Engineering Equipment](#), [Strength Of Material Lab](#)

Product Page:

<https://www.labappara.com/product/rockwell-brinel-hardness-tester/>

Product Description

Rockwell & Brinel Hardness Tester

Metal testing requires various stages and Rockwell Cum Brinel Hardness Tester is essential for calculating the hardness of metal. The process is automated and you can rely on the automatic load selection in accordance to the varying sample. They can measure the hardness of metal and its alloys irrespective of their natural hardness and softness.

Features:

- Multiple scales
- Varying loads
- Automated process
- Resistant to abrasion.
- Robust construction.
- Smooth finish.

Specifications:

- Scales- a,b,c & brinell
- With brinell microscope and all essential accessories

Automatic load selection with automatic zero setting dial gauge

Rockwell 60, 100 and 150 kgf

Brinell 187.5 kgf

With ball indenter