

Rota Rod



Product Categories: [Pharmaceutical Equipments](#)

Product Tags: [Laboratory Equipment](#)

Product Page: <https://www.labappara.com/product/rota-rod/>

Product Description

Rota Rod - Pharmaceutical Equipment

Rota Road Apparatus - Latest with Digital Counter for determining neuro-toxicity, muscle tone, balance and motor co-ordination in rats and mice. Consisting of two or three or four compartments of 75 mm width each with a rotating rod of 25 mm diameter having speed of 5,10,15,20 and 25 rpm with time interval counters in each compartment.

The Rota rod performance test is a performance test based on a rotating rod with forced motor activity being applied, usually by a rodent.

The test measures parameters such as riding time (seconds) or endurance

Two Compartment

Three Compartment

Four Compartment

Specification:

Sturdy M.S. frame covered with powder coating gives a fine look to apparatus.

Digital timer control panel.

Driven by a D.C. drive in variable speed to the apparatus.

The compartment is divided by a transparent acrylic sheet.

Features:

Best grade materials

Compatible

High performance