

Round Test Tube Stand



Product Categories: [Clamps, Holders Racks Stands](#), [Laboratory Plastic Ware](#)

Product Tags: [Laboratory Equipments](#), [Laboratory plastic ware](#)

Product Page: <https://www.labappara.com/product/round-test-tube-stand/>

Product Description

Round Test Tube Stand - Scientific Laboratory Plastic ware

Test tube stand is a laboratory equipment that is used to hold upright multiple test tubes at the same time. Our Round Test Tube Stand is Space Saver.

A handy, little, space saver, this circular twelve place Polypropylene Test Tube Stand has four 25 mm dia. Holes & eight 19 mm dia. Holes in its top plate. The base has matching hemispherical wells & 16 vertical pins for drying tubes surrounded by a trough to catch drippings. This autoclavable rack can be easily disassembled for cleaning purpose.

Material : Polypropylene

Particulars	Packing
19 mm & 25 mm	2 Pcs