

SCR Characteristics Apparatus



Product Categories: [Characteristics And Application](#), [Electronics](#)

Product Page: <https://www.labappara.com/product/scr-characteristics-apparatus/>

Product Description

SCR Characteristics Apparatus

Features :

Instrument comprises of two DC Regulated Power Supplies 0-30VDC/150mA & 0-1VDC/150mA, three round meters for voltage & current measurement, SCR 2P4M mounted behind the panel, connections of Supplies, Meters & SCR brought out at 4mm Sockets.