


## Semi Fowler Bed



**Product Categories:** [Hospital Equipment & Furniture](#)

**Product Tags:** [hospital](#), [medical](#)

**Product Page:** <https://www.labappara.com/product/semi-flower-bed/>

### Product Description

Semi Fowler Bed - Hospital Equipment & Furniture

Mechanically operated two sectional, uniformly perforated CRCA sheet top beds.

ABS head and foot bows / panels

Strong "H" Shaped legs on rubber stumps.

One Fix lever handles at foot end to obtain head raise upto 80deg.

**Overall Size** - 208L X 92W X 60H