

## Semi Micro Burner



**Product Categories:** [Chemistry Equipments](#), [Lab Supplies](#)

**Product Tags:** [chemistry lab instruments](#), [science lab equipment](#)

**Product Page:** <https://www.labappara.com/product/semi-micro-burner/>

### Product Description

Semi Micro Burner- Laboratory Equipment

Laboratory Equipment that produces a single open gas flame.

It is used for heating, sterilization, and combustion.

The gas can be natural gas (which is mainly methane) or a liquefied petroleum gas, such as propane, butane or a mixture of both.

### Feature

Wide & stable base provide.

Brass parts in iron base.

For L.P.G of petrol.

Nicely Packed in box