

## Senior Microscope 6"



**Product Categories:** [Audio- Visuals](#), [Educational equipments](#)

**Product Tags:** [biology lab equipment](#), [chemistry lab instruments](#)

**Product Page:** <https://www.labappara.com/product/senior-microscope-6/>

### Product Description

#### Senior Microscope 6" - **Scientific Laboratory Educational Equipment**

The senior microscope is defined as the shortest distance between two points on a specimen that can still be distinguished by the observer or camera system as separate entities.

#### **Feature -**

Two objective 10X, 45X.

100 mm Dia screen Projection Microscope.

One eye piece 10X.

Fitted with 12v- 21w halogen lamp.

Mechanical Stage.

With intensity control knob.

Nicely packed in wooden box.

Microscope Body is heavy duty, inclined monocular type with sturdy base.

Equipped with Coarse & Fine Focusing systems with Quadruple revolving Nose Piece.

Plain stage (120x120mm) with graduated Mechanical Stage for X-Y movements of specimen slide.

Built in light by 6V-20W halogen lamp with continuously variable light intensity control.

The whole system of condenser & light arrangement is moveable up & down by rack & pinion to get the optimum results.

Projection Doom of screen size 150 mm is rotatable through 360°.

Objectives : Achromatic 10x & 20x/40x (SL)

Projection Eye Piece : 10x (Huy)

Packed in a thick and attractive plywood box.