

Senior Projection Microscope



Product Categories: [Audio- Visuals](#), [Educational equipments](#)

Product Tags: [chemistry lab instruments](#), [physics lab instruments](#)

Product Page:

<https://www.labappara.com/product/senior-projection-microscope-2/>

Product Description

Senior Projection Microscope - Scientific Laboratory Educational Equipment

The senior projection microscope is defined as the shortest distance between two points on a specimen that can still be distinguished by the observer or camera system as separate entities. A test improved model for accurate measurement of micro structures, solid state light control with digital intensity meter and special pointer system . Having 200 mm graduated screen rotatable 360 with graduation around the screen and vernier reading to 0.1 degree. Equipped with low Halogen lamp and imported Fresnel lens behind the screen for uniform projection throughout the screen with achromatic objectives as follows:

With Indian Optics

With optics having imported lenses. Extra charges for stage

Objectives 5 x, 10 x, 20 x, 40 x, 60 x, 5 Nos.

Projection Range 100 x to 1200 x

Voltage 100 watts

Feature -

Compact dust proof ,construction projection

Screen of 150 mm dia.

With provision of wall projection,

12 v - 100 w halogen lamp.

Equipped with loop and arrow pointer.

Supplied with mechanical stage and two objectives 10 x 740 x.

Eye piece packed in card board box.