

Senior Projection Microscope



Product Categories: [Microscope](#), [Microscopes](#)

Product Tags: [biology lab equipment](#)

Product Page: <https://www.labappara.com/product/senior-projection-microscope/>

Product Description

Senior Projection Microscope - Laboratory Equipment

It is an advanced versatile projection microscope for measurement of specimen, equipped with bright Halogen Lamp 12 V, 100 Watts solid state variable light control arrangement workable on 220V A.C mains, inter locked with electronic fuse system. Quadruple revolving nose piece, approximate 200 mm aid.

Magnification from 125X to 1000X

Objectives. 5X,10X,20X & 40X