

Senior Student Microscope



Product Categories: [Microscope](#), [Microscopes](#)

Product Tags: [biology lab equipment](#)

Product Page: <https://www.labappara.com/product/senior-student-microscope/>

Product Description

Senior Student Microscope - Scientific Equipment

Olympus (HSA) type eyepiece 10 x, 15 x and objectives 10 x & 40 x, having condensers packed in box. Focusing rack and pinion knob provided for sensitive focusing.

Knobs separate coarse and fine motion knobs with very precise & accurate graduated fine motion reading to 0.002 mm, for fine adjustment and a stopper to prevent damage to glass micro slides and objective lenses.

Illumination standard plano concave mirror 50 mm dia. In fork for natural light, replaceable by a sub-stage microscope lamp.