

Separatory funnel



Product Categories: [Funnels](#), [Laboratory Plastic Ware](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#), [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/seperatory-funnel/>

Product Description

Separatory funnel - Scientific Laboratory Equipment

Separatory Funnel is made of Polypropylene, giving it excellent contact clarity and strength. This pear shaped, autoclavable separating funnel has a leakproof threaded screw cap. The Polypropylene stopcock, fitted with a self lubricating PTFE plug, is push fitted to the stem of the separatory funnel. A separating funnel takes the shape of a cone with a hemispherical end. It has a stopper at the top and stopcock (tap), at the bottom.