

## Single Beam UV-VIS Spectrophotometer



**Product Categories:** [Analytical Lab Instruments](#)

**Product Tags:** [Analytical Laboratory Equipment](#)

**Product Page:**

<https://www.labappara.com/product/single-beam-uv-vis-spectrophotometer/>

### Product Description

Single Beam UV-VIS Spectrophotometer - Analytical Laboratory Equipment

Range 200-1000 nm with facility for Automatic single beam uv-vis spectrophotometer.

Concentration, % Transmission and absorbency Measurement.

Wavelength resolution 5nm 3 and half digit LCD display With cuvettes a set of 1 pair.

Cover 195 to 1000 nm, 4 Position Cuvette Holder RS232C Computer Interface.

Highly Accurate & Stable High wavelengths Precision Compact & Good Looking Shape.

A versatile, Rugged, Easy to operate and Accurate Instrument to determine

## Absorbance, % Transmission and Concentration of Liquids.