

Single Nose Microscope



Product Categories: [Microscope](#), [Microscopes](#)

Product Tags: [biology lab equipment](#)

Product Page: <https://www.labappara.com/product/single-nose-microscope-2/>

Product Description

Single Nose Microscope - Laboratory Equipment

With one eyepiece 10 x and one objective 10 x

Total Magnification: 100 x Plano concave mirror on gamble mount for illumination, Gives Clear Images, Easy To Operate, Ideal for School & Kids. Well packed in thermometer box.

Description:

Sand casted body Inclination 90° on Horse shoe base. Focused by rack and pinion. Total Height 280 mm & Mechanical Tube Length 160 mm. Fixed square stage size 100 mm x 100 mm with stage clips. Nicely packed in Thermocol Box for packing eye piece.

Features:

Smooth functionality

Long working life

Minimum maintenance

Optical Combination:

Eyepiece 10 x

objective 10 x

Magnification: 100 x