

SINGLE NOSE MICROSCOPE



Product Categories: [Biology Lab](#), [Educational equipments](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#), [Manufacturer of Laboratory equipment](#), [physics lab instruments](#), [science lab equipment](#)

Product Page: <https://www.labappara.com/product/single-nose-microscope/>

Product Description

SINGLE NOSE MICROSCOPE - Biology Laboratory Equipment

Sand casted body Inclination 90° on Horse shoe base. Focused by rack and pinion. Total Height 280 mm & Mechanical Tube Length 160 mm. Fixed square stage size 100 mm x 100 mm with stage clips. Nicely packed in Thermocol Box for packing eye piece.

Features:

Smooth functionality

Long working life

Minimum maintenance