

## Sintered Wire



**Product Categories:** [Laboratory Glassware](#), [Sintered Ware Borosilicate Glass](#)

**Product Tags:** [glassware](#)

**Product Page:** <https://www.labappara.com/product/sintered-wire/>

### Product Description

Sintered Ware - Laboratory Glassware Equipment

Sintered ware Tube with Sintered Disc dia. Sealed in Center, Porosity G-0 or G-1 or G-2 or G-3 G-4

Sintered ware is used for the filtration of liquids and gases in the laboratory. It incorporates a porous glass disc as a filter media, which is non corrosive and reusable. It is also used for gas washing, dispersion and absorption. Sintered disc is manufactured by crushing brand glass, powdering, cleaning, separating into various mesh sizes and then fusing together in the form of a disc.

Available in Diameter - 10, 20, 30, 45, 50, 90