

## Slide Storage Rack



**Product Categories:** [Laboratory Plastic Ware](#), [Slide Storage](#)

**Product Tags:** [biology lab equipment](#), [chemistry lab instruments](#), [Laboratory Equipment](#), [Laboratory plastic ware](#)

**Product Page:** <https://www.labappara.com/product/slide-storage-rack/>

### Product Description

Slide Storage Rack - Scientific Laboratory Plastic ware

Slide Storage Rack can hold 100 slides which are recognizable through the numeric indexing. Adding to the existing range of slide accessories, this slide storage rack made of ABS can hold 100 slides, recognizable through the numeric indexing. Spacing between the slots provide adequate backward & forward movement of slides for convenient recognition. The transparent cover provide excellent view of the slides stored inside. It is widely used for storing slides to prepare samples in laboratories.

#### **Features:**

Crack resistance

Light weight

Smooth finish

Durability

Slide Storage Rack - 100 Slides