

Slide Warming Table



Product Categories: [General Chemistry Labware](#), [Hot Plate](#), [Lab Heating-Equipment](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/slide-warming-table/>

Product Description

Slide Warming Table - Scientific Equipment

Thermostatically controlled uniform temp. From room to 800 C top made of anodized aluminum / stainless steel.

Slide Warming Table is very useful for heating slides to a uniform constant temperature below the melting point of wax

Top surface made of Stainless Steel

The surface temperature is controlled from ambient to 70°C

Supplied complete with pilot lamp, cord and plug but without thermometer

To work on 220/230 Volts A.C. Supply

The body is made out of thick mild steel duly finished in white stoving enamel/powder coated paint with mat finished color combinations and top is made of highly polished aluminum sheet.

Top Plate Sizes Length x Width

30 x25 cm

12" x 10"

40 x 25 cm

16" x 10"

30 x 45 cm

18" x 12"

60 x 45cm

24" 18"