


Small Animal Cage



Product Categories: [Laboratory Equipments](#), [Pharmaceutical Equipments](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/small-animal-cage/>

Product Description

Small Animal Cage - Laboratory Equipment

A rectangular plastic container with a stainless steel grilled top for ventilation, suitable for the species under study. The grill is so designed for placing the water bottle for feeding. Available in 2 different cage and feeding bottle sizes .

Poly Carb, Bod Clear - 370 cc

Dimensions: 267 x 207 x 140 mm approx

Floor Area: 370 cm²

Auto-clavable at 1210 C

Able to house 6 small mice (<20 g) or up to 3 large mice (>30 g) as per EU Guidelines

Made of steel body duly powder coated or stainless steel.

With feeding bottle and feeding bowl.

The removable waste collecting tray.