

SOLAR FAN (100mm)



Product Categories: [Educational equipments](#), [Geography & Working Models](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/solar-fan-100mm/>

Product Description

SOLAR FAN (100mm) - Scientific Laboratory Educational Equipment

When solar cell is placed in sunlight it converts solar energy into electrical energy, which in turn energizes the motor and thus conversion of solar energy into electromechanical energy.

This is a great introduction into the workings and benefits of solar energy and it is a fun and educational project for the classroom.