

Solar Meter



Product Categories: [Educational equipments](#), [Geography & Working Models](#)

Product Tags: [Eduactional Equipments](#), [Lab Equipments](#)

Product Page: <https://www.labappara.com/product/solar-meter/>

Product Description

Solar Meter - Scientific Laboratory Educational Equipment

Instrument is used to measure the power density of sunlight. Moving coil DC meter housed in non-breakable ABS plastic case with calibrated scale expressed in mW/cm^2 to show full sun intensity level and encapsulated with Solar cell.

With the help of our professionals, LABAPPARA have been able to provide an excellent array of Solar Meter which is sourced from the leading vendors.

Our customers can avail from us a complete range of Solar Meter.

It helps to demonstrate, analyze and calculate the solar energy.

Features:

Excellent service back up

Zero maintenance

Durable