

## Sonometer (Aluminum)



**Product Categories:** [Physics Equipments](#), [Sound And Waves Oscillations](#)

**Product Tags:** [physics lab instruments](#)

**Product Page:** <https://www.labappara.com/product/sonometer-aluminum/>

### Product Description

#### **Sonometer (Aluminum)**

A Sonometer is a device for demonstrating the relationship between the frequency of the sound produced by a plucked string, and the tension, length and mass per unit length of the string. It is used in School, College, University, Research Laboratories.

#### **Feature:**

Sonometer Aluminum Soft wood nicely polished with metal fitting.

Sonometer Aluminum With wooden riders.

And provide with 3 different wire 100 cm long.

Sonometer Aluminum fitting is without wt.